Chapter ATCP 40

APPENDIX B

HUMIC SUBSTANCE EXTRACTION TEST

The analytical method for determining the percentage content of a humic substance in a product shall follow a colorimetric test of a 0.5N NaOH extraction from the formulated product. The methodology is as follows:

Scope

This method may be used for solid and liquid samples containing 0.5% or more humic acid. The humic acids are dissolved by treatment with 1N sodium hydroxide and then precipitated with hydrochloric acid.

Equipment and Reagents

- 1. Centrifuge
- 2. 100 mL Screw Cap Centrifuge Bottles
- 3. 1 N NaOH
- 4. 1% NaOH (10 g/liter)
- 5. Conc. HCl
- 6. 100°C Drying Oven
- 7. Mechanical Shaker

Analysis

- 1. Weigh appropriate size sample into a 100 mL wide mouth screw top bottle to give close to 0.2 gms of dry humic acid ppt.
- 2. Add 50 mLs 1 N NaOH, seal tightly.
- 3. Shake on mechanical shaker for 1.5 hours for solids, 30 minutes for liquids.
- 4. Rinse the cap with 5 mLs 1% NaOH.

- 5. Centrifuge for 25 minutes at 2000 rpm.
- 6. Decant supernatant liquid into a second weighed bottle.

ATCP 40 Appendix B

- 7. Add 10 mLs 1% NaOH to first bottle, shake vigorously, centrifuge again.
- 8. Add the supernatant liquid to the second centrifuge bottle.
- 9. To the combined extracts in the second bottle, add conc. HCl until the pH is adjusted to between 1 and 2.
- 10. Centrifuge the sample for 25 minutes at 2000 rpm.
- 11. Carefully decant the liquid and discard.
- 12. Add 25 mLs distilled water (previously adjusted to pH 1-2 with HCl) to the bottle with ppt., shake vigorously to free all ppt. from bottom and centrifuge again.
- 13. Again, carefully decant the liquid and discard.
- 14. Repeat steps 12 and 13 two more times.
- 15. Dry the bottle with humic acid overnight at 100° 110° C.
- 16. Cool in dessicator and weigh.

Calculations

% Humic Acid = Weight dried residue X 100
Sample weight

References

John Husler, University of New Mexico, Department of Geology, Albuquerque, New Mexico

C.A. Black, <u>Methods of Soil Analysis</u> Part 2, American Society of Agronomy, Inc., Madison, Wisconsin, 1965.